#### REPORT 1

- 2
- 3

# Peripherally-derived LGI1-reactive monoclonal antibodies cause epileptic seizures in vivo

Manoi Upadhya,<sup>1</sup> Toni Kirmann,<sup>2</sup> Max A. Wilson,<sup>1</sup> Christian M. Simon,<sup>2</sup> Divya Dhangar,<sup>1</sup> 4

Christian Geis,<sup>3</sup> Robyn Williams,<sup>4,5</sup> Gavin Woodhall,<sup>1</sup> Stefan Hallermann,<sup>2</sup> Sarosh R. Irani<sup>4,5</sup> 5

6

and Sukhvir K. Wright<sup>1,6</sup>

#### Abstract 7

One striking clinical hallmark in patients with autoantibodies to leucine-rich glioma 8 inactivated 1 (LGI1) is the very frequent focal seizure semiologies, including faciobrachial 9 dystonic seizures (FBDS), in addition to the amnesia. Polyclonal serum IgGs have 10 successfully modelled the cognitive changes in vivo but not seizures. Hence, it remains 11 unclear whether LGI1-autoantibodies are sufficient to cause seizures. 12

- We tested this with the molecularly precise monoclonal antibodies directed against LGI1 13 (LGI1-mAbs), derived from patient circulating B cells. These were directed towards both 14 major domains of LGI1, LRR and EPTP and infused intracerebroventricularly over 7 days 15 into juvenile male Wistar rats using osmotic pumps. Continuous wireless EEG was recorded 16 from a depth electrode placed in hippocampal CA3 plus behavioural tests for memory and 17 hyperexcitability were performed. Following infusion completion (Day 9), post-mortem brain 18 slices were studied for antibody binding and effects on Kv1.1. 19
- The LGI1-mAbs bound most strongly in the hippocampal CA3 region and induced a 20 significant reduction in Kv1.1 cluster number in this subfield. By comparison to control-Ab 21 injected rats video-EEG analysis over 9 days revealed convulsive and non-convulsive seizure 22 activity in rats infused with LGI1-mAbs, with a significant number of ictal events. Memory 23 was not impaired in the novel object recognition test. 24
- 25 Peripherally-derived human LGI1-mAbs infused into rodent CSF provide strong evidence of direct in vivo epileptogenesis with molecular correlations. These findings fulfill criteria for 26 LGI1-antibodies in seizure causation. 27
- 28

#### Author affiliations: 29

1 Institute of Health and Neurodevelopment, School of Health and Life Sciences, Aston 30

University, Birmingham, B4 7ET, UK 31

> © The Author(s) 2024. Published by Oxford University Press on behalf of the Guarantors of Brain. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (https://creativecommons.org/licenses/by/4.0/), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.

- 1 2 Faculty of Medicine, Carl-Ludwig-Institute of Physiology, Leipzig University, Leipzig
- 2 04103, Germany
- 3 3 Department of Neurology, Section Translational Neuroimmunology, Jena University
- 4 Hospital, Jena 07747, Germany
- 5 4 Oxford Autoimmune Neurology Group, Nuffield Department of Clinical Neurosciences,
- 6 University of Oxford, OX3 9DU, Oxford, UK
- 7 5 Departments of Neurology and Neuroscience, Mayo Clinic, Jacksonville, Florida 32224,
- 8 USA
- 9 6 Birmingham Women's and Children's Hospital NHS Trust, Birmingham, B4 6NH, UK
- 10
- 11 Correspondence to: Dr Sukhvir Wright
- 12 Institute of Health and Neurodevelopment, School of Health and Life Sciences, Aston
- 13 University, Birmingham, B4 7ET, UK
- 14 E-mail: <u>s.wright5@aston.ac.uk</u>
- 15
- 16 **Running title**: LGI1 monoclonal antibody epileptogenicity
- 17

18 Keywords: LGI1-Ab encephalitis; faciobrachial dystonic seizures (FBDS); autoimmune19 associated epilepsy; seizures; Kv1.1

20

Abbreviations: AMPAR α-amino-3-hydroxy-5-methyl-4-isoxazole propionic acid CSF
 cerebrospinal fluid EPTP epitempin repeat FBDS Faciobrachial dystonic seizures IgG
 Immunoglobulin LE Limbic encephalitis LGI1 leucine-rich glioma inactivated 1 LRR
 leucine-rich repeat NOR Novel object recognition PSBB Post seizure behavioral battery ROI
 region of interest

26

# 27 Introduction

Autoantibodies to leucine-rich glioma inactivated 1 (LGI1) are identified mainly in older males with a variety of distinctive and frequent focal seizure semiologies, including faciobrachial dystonic seizures (FBDS), piloerection and thermal seizures. In addition, many of these patients develop profound amnesia<sup>1,2</sup>. FBDS are characterised by brief and numerous seizures characterised predominantly by contractions of the arm and ipsilateral face<sup>3</sup>. If unrecognised and untreated, FBDS can evolve into a limbic encephalitis (LE) with temporal lobe and tonic-clonic seizures<sup>3</sup>,<sup>4</sup>, whereas, early treatment of FBDS with immunotherapies
 typically dramatically reduces seizure frequencies and can prevent the development of LE,
 thereby avoiding additional disability<sup>4</sup>.

Medial temporal lobe changes on MRI are a common finding in patients with LGI1antibodies, and the subsequent atrophy predominantly affects the CA3 region<sup>5</sup>. Ictal EEG
changes have been inconsistently reported, however, the most commonly recognised patterns
are generalised attenuation, electrodecremental changes <sup>2, 6, 7</sup>, simultaneous frontal EEG and
EMG unilateral slow and infraslow wave preceding the contralateral tonic-dystonic seizure<sup>8</sup>,
as well as frequent *subclinical* temporal lobe seizures <sup>2, 8, 9</sup>.

Polyclonal patient serum-derived IgGs and LGI1-reactive monoclonal antibodies cause 10 neuronal hyperexcitability and even epileptiform activity in vitro, with AMPAR and/or Kv1.1 11 downregulation postulated as the most likely underlying molecular mechanisms<sup>10,11</sup>. 12 However, despite seizures being a key clinical hallmark of this disease, animal models using 13 LGI1-antibodies have failed to recapitulate this pathognomonic feature in vivo <sup>12,13</sup>. This 14 observation calls into question the direct epileptogenicity of LGI1-antibodies, particularly as 15 the use of polyclonal human serum means that non-LGI1 reactivities may mediate observed 16 outcomes. 17

Here, we use human-derived monoclonal antibodies with specific reactivities to either the leucine-rich repeat (LRR) or the epitempin repeat (EPTP) domain of LGI1<sup>14</sup>, to produce a passive transfer rodent model with spontaneous epileptic seizures and video-EEG recorded ictal changes.

22

# 23 Materials and methods

# 24 Antibody preparation

Two LGI1-specific mAbs, targeting either the LRR (mAb2) or EPTP (mAb13) domains were
 produced, as described previously <sup>14</sup>.

## 27 In vivo experiments

## 28 Animals

Sixteen postnatal day 21 (P21) male Wistar rats, weighing 52-60 g were used for experiments
described. Animals were housed in temperature- and humidity-controlled conditions with a
12 h/12 h light/dark cycle and allowed free access to food and water. All procedures were

compliant with current UK Home Office guidelines as required by the Home Office Animals
 (Scientific Procedures) Act 1986 and carried out under the authority and procedural approval
 of a UK Home Office approved project license and in line with ARRIVE guidelines. The
 Aston Bioethics Committee, University of Aston, Birmingham, UK granted local ethical
 approval for the study. Full experimental timeline shown in Figure 1A.

# 6 EEG implantation and osmotic pump surgery

The rats were implanted, as previously described, with subcutaneous transmitters for in vivo 7 EEG recordings (A3028B-DD subcutaneous transmitters, 90-mm leads, OpenSource 8 Instruments (OSI)) via unilateral depth electrode (W-Electrode (SCE-W), OSI) in left 9 hippocampus (CA3, 3.5 mm lateral, 3.8 mm caudal, depth 2.3 mm), and subcutaneous 10 11 osmotic pumps (model 1007D, Azlet) (volume 100 µl, flow rate 0.5 µl/hr, duration 7 days, primed overnight) attached to bilateral cannulae (328OPD-3.0/SPC; PlasticsOne), implanted 12 into the lateral ventricles (±1.5 mm lateral and 0.6 mm caudal) <sup>15-17</sup>. A reference EEG 13 electrode was implanted on the contralateral skull surface (3.5 mm lateral, 3.8 mm caudal); 14 the cannula and skull electrodes were secured with dental cement <sup>15,16</sup>. 15

## 16 EEG data analysis

EEG data (wide band pass 0.2-160 Hz sampled at 512 samples per second) was collected and 17 recorded using Neuroarchiver software (OSI) from wireless transmitters implanted in freely 18 moving animals placed in a custom-built Faraday cage with aerial (OSI) as previously 19 described <sup>15-17</sup>. For automated ictal event detection, video-EEG matching was used to identify 20 ictal EEG events. The Event Classifier (OSI) was used to classify segments (1 second) of 21 22 EEG according to program metrics (power, coastline, intermittency, coherence, asymmetry, spikiness) creating clusters of similar events when plotted. This generated a library of ictal 23 events that allowed fast identification of abnormal EEG events by automated comparison 24 (http://www.opensourceinstruments.com/Electronics/A3018/Seizure\_Detection.html). 25

26 Powerband analysis was carried out using a custom-designed macro.

# 27 Behavioural testing

## 28 Novel Object Recognition (NOR)

The novel object recognition test is used to test cognitive performance in autoimmune encephalitis models <sup>13,14</sup>. Animals were habituated to the test arena 24 hours before training and tests. The time interacting with the objects was measured and the NOR index (time interacted with novel object/ total time of interaction with both the objects) calculated. Three tests were carried out at regular intervals during the 9-day recording period. All the observations were carried out by the experimenter blind to the treatments and analysed using the Ethovision software. The frequency of rats to approach the object, time spent with the objects was measured for NOR index, and the distance travelled, and velocity was measured as a measure of anxiety score.

## 7 Post-Seizure Behavioural Battery (PSBB)

8 The Post-Seizure Behavioural Battery (PSBB) was performed as previously described to 9 assess hyperexcitability and other behavioural indices suggestive of epilepsy <sup>16,18</sup>. Two simple 10 and non-stressful tasks, "touch" and "pick-up" tasks constituted the PSBB and were 11 performed at three time-points. PSBB scores were calculated by taking the product of the task 12 scores ('touch x pickup').

13

## 14 Local field potential recordings

### 15 **Local field potential (LFP) recordings:**

LFP recordings were performed and analysed as previously described.<sup>15,16</sup> Briefly, on day 9 16 17 immediately after the behavioral experiments, rats were anaesthetized using isoflurane and 18 following the loss of consciousness, pentobarbital (60 mg/kg, SC) and xylazine (10 mg/kg, IM) injected. Transcardial perfusion was then performed using ice-cold artificial 19 cerebrospinal fluid (aCSF). Animals were decapitated, and the brain was removed. Brain 20 slices were prepared using a vibratome (Campden Instruments, UK) at 450µm for LFP 21 recording and 350µm for fluorescence intensity measurements. LFP recordings were 22 assessed for spike activity using Spike2 software (CED). Root mean square (RMS) amplitude 23 of each recording was calculated, and events with an amplitude greater than five-fold the 24 RMS amplitude were considered a "spike". 25

26

# 27 Immunofluorescence, immunohistochemistry and image analysis

- 28 The immunofluorescence study was performed as previously reported  $^{15,16}$ .
- 29 Further detailed methods available in Supplementary information.
- 30

## 1 Statistical analysis

For behavioural experiments repeated-measures 2-way ANOVA with Bonferroni multiple
comparison test and non-parametric Mann-Whitney T-test were used to assess differences.
For in vivo experiments and immunohistochemistry, two-tailed Mann-Whitney T-test were
used. P<0.05 was considered as significant. Statistical analysis was conducted using</li>
GraphPad Prism 8 (GraphPad Software Inc).

7

# 8 **Results**

# 9 LGI1 antibodies bind to rat hippocampus and cause reduction of

## 10 Kv1.1 clusters

Seven-day osmotic pumps delivered the LGI1-mAbs (n=10; five animals with LRR and five
with EPTP specific mAbs) and control Abs (n=6) into the lateral cerebral ventricles of
P21Wistar rats.

After interventricular infusion, both LGI1-mAbs (LRR and EPTP specific) bound to rodent 14 hippocampus (Figure 1B), confirming previous studies in mice <sup>14</sup>. By comparison to control, 15 non-brain-reactive Abs, the fluorescence intensities of bound LGI1-mAbs were highest in the 16 CA3 and dentate gyrus (Figure 1C; Supplementary Figure 1A) This increased binding was 17 18 also seen in the contralateral hemisphere (without depth electrode; Supplementary Figure 1B). As demonstrated with prior use of these antibodies, the LRR mAbs were strongly 19 retained in the infused hippocampus, EPTP showed slightly less residual binding overall on 20 both sides (Supplementary Figure 1 A, B).<sup>14</sup> In the CA3 region, super-resolution light 21 microscopy using stimulated emission depletion (STED) of 20-µm-thick cryosections 22 revealed that the number of Kv1.1 clusters was significantly lower in the LGI1-mAb infused 23 group (pooled analysis of LRR and EPTP) as compared to controls (Figure 1D,E), and non-24 25 significantly trended to those which co-localised with the glutamatergic synapse marker vGlut1 (Figure 1F,G). 26

27

# LGI1-mAbs cause clinical and subclinical spontaneous epileptic seizures *in vivo*

To investigate potential epileptogenic effects, brain activity was recorded continuously during
 LGI1-mAb infusion using wireless EEG transmitters from a depth electrode placed in the

1 CA3 region of hippocampus. One rodent infused with EPTP specific LGI1-mAb was culled 2 early due to excess lacrimation, and in one rodent infused with LRR specific mAb, the EEG 3 recording was suboptimal for analysis but the animal data were still included for behavioural 4 testing and post-mortem immunohistochemistry. There was no mortality or other significant 5 morbidity following seizures and documented welfare scores remained low throughout the experimental period indicating minimal pain or distress. There was ictal activity recorded in 6 one control animal within 12 hours of surgery, most likely due to acute effects of 7 surgery/anaesthesia. The PSBB score remained <2 for this animal throughout the 8 experimental period with no further ictal activity noted on EEG. In the LGI1-mAb infused 9 rats, the frequency, duration and recurrence of ictal events as identified through automated 10 seizure detection <sup>15-17</sup> in each animal treated is shown in Figure 2E,F. The recorded epileptic 11 and behavioural activity included convulsive and non-convulsive ictal events (Figure 2A, E, 12 F; Suppl. videos 1,2,3). We observed a significantly higher total number of ictal events in the 13 LGI1-mAb infused animals compared to controls during the 9 days of EEG recording period 14 (Figure 2B). EEG coastline was higher in LGI1-mAb infused animals, due to higher 15 16 prevalence of relatively high-amplitude epileptiform events (Figure 2C). This cohort also showed higher power in all EEG frequency bands. Similar spectral power changes compared 17 to controls are seen in other autoimmune-encephalitis-associated seizure models and models 18 19 of status epilepticus, reflecting the presence of recurrent spontaneous convulsive and nonconvulsive seizures <sup>19-22</sup>(Figure 2D). Following completion of the antibody infusion, LGI1-20 mAb and control-Ab infused animals were terminally anaesthetized and brains extracted for 21 22 acute in vitro brain slice recording. Local field potentials were recorded in CA3 and CA1 of the hippocampus, spike activity was not significantly different from control slices 23 (Supplementary Figure 1C). This is concordant with the minimal in vivo ictal activity 24 recorded by Day 9 (Figure 2E, F). 25 26

# LGI1-mAbs cause enhanced startle response and aggressive behavior without effect on cognition

The post-seizure behavioural battery (PSBB) was used to monitor behavioural changes (enhanced startle response and aggression) consistent with the development of spontaneous recurrent seizures <sup>16,18</sup>. Animals were tested three times over 9 days using touch and pick-up tests, with a score > 10 indicating significant behavioural change. The LGI1-mAb infused animals scored significantly highly at all three time-points compared to the control-Ab infused animals where the scores remained consistently low (Figure 3A; Supplementary
Figure 1). There were no differences in the NOR index (calculated by dividing the time spent
with novel object to the entire duration of time spent with both the objects), locomotion, and
distance travelled and velocity, during NOR across three timepoints (day 1, day 5 and day 9;
Figure 3B-D; Supplementary Figure 1).

6

# 7 **Discussion**

LGI1-antibody encephalitis patients exhibit focal and tonic-clonic seizures as well as 8 cognitive impairment. Published animal models were based on application of serum-IgGs to 9 successfully model the cognitive changes in vivo, but seizures were not seen <sup>12,13</sup>. Here, we 10 used highly specific patient-derived LGI1-mAbs, which lack the potential additional 11 reactivities potentially present in serum. After infusion into the CSF of juvenile Wistar rats, 12 these LGI1-mAbs led to seizures in vivo. Similar to patient phenotypes, where epileptiform 13 discharges are the most frequent EEG finding<sup>19</sup>, we recorded clinical and subclinical 14 seizures, with ictal EEG changes. These findings support the direct epileptogenicity of LGI1 15 antibodies and allow them to fulfill Wietebsky's postulates of pathogenicity <sup>20</sup>. With 16 17 continuing controversy over whether LGI-Ab associated seizures are in fact a paroxysmal movement disorder and not ictal<sup>21</sup>, this study confirms the causal link between the molecular 18 interactions affecting Kv1.1 mediated by LGI1-Abs and seizure activity in vivo. 19

LGI1 is a secreted neuronal protein that interacts with the ADAM22 and ADAM23 20 transmembrane metalloproteases forming a trans-synaptic protein complex that facilitates 21 excitatory synaptic transmission<sup>22</sup>. LGI1 consists of an N-terminal LRR domain and a C-22 23 terminal EPTP domain. The LRR-directed monoclonal antibody used in this study binds the ADAM22/23-docked LGI1 complex in vitro, causing internalisation and hence disruption of 24 the membrane stabilising trans-synaptic complex<sup>14</sup>. This pathomechanism is similar to that 25 demonstrated by a CSF-derived LGI1-specific mAb that also caused epileptiform changes in 26 vitro<sup>10</sup>. In vitro primary neuronal cultures also incubated for 7 days with LRR and EPTP 27 specific mAbs showed that LRR-directed mAbs could directly affect neuronal excitability, 28 but this effect was less pronounced in EPTP-specific mAb exposed neurons<sup>23</sup>. In contrast to 29 the action of LRR mAbs in vitro, the EPTP mAbs are reported to exert their effects by 30 *directly* inhibiting the docking of LGI1 to the ADAM proteins <sup>14</sup>. In our study, EEG 31 recordings have demonstrated increases in network and neuronal excitability in hippocampal 32 CA3 in vivo by both LRR and EPTP domain directed LGI1 autoantibodies <sup>10, 23</sup>. Studies have 33

shown that LRR- and EPTP-specific antibodies co-occur in patients with LGI1-antibody encephalitis but the numbers of recorded rats with either mAb subclass do not allow for statistical comparison of epitope-specific clinical features (Supplementary Figure 1) <sup>14</sup>. The goal of future studies is to analyse our *in vivo* EEG using computational modelling to evaluate different hypotheses pertaining to network hyperexcitability and seizures by action of epitope specific LGI1-antibodies and to differentiate epitope-specific effects <sup>15, 24</sup>.

EPTP specific mAbs are reported to have no effect on cognitive performance when injected 7 into hippocampi of mice <sup>14</sup>: this may explain why our pooled NOR results from both LRR 8 and EPTP specific mAbs infused rats were unremarkable. Additionally, in contrast to 9 previous in vivo passive transfer models, our antibody infusion time was much shorter: 7 10 days, compared to 14 days. Patients with LGI1-Ab encephalitis frequently develop seizures 11 before the onset of memory disturbance<sup>25</sup> and this could explain our behavioural findings in 12 this rodent model. Lastly, the genetic background of rodents can impact on seizure 13 susceptibility <sup>12, 17, 26</sup>; here juvenile Wistar rats were used as they have proven to be effective 14 in modelling autoimmune-associated seizures and epilepsy in vivo <sup>15, 16</sup>. 15

16 The goal of this study was to prove epileptogenicity of LGI1 mAbs in vivo using mAbs with known pathogenic effects in vitro. We also took this opportunity to generate some insights 17 into identifying a pre- or post- synaptic mechanism for seizure generation. Previous longer-18 19 term infusion (14 days) in vivo studies using patient-derived LGI1-mAbs that target both the LRR- and EPTP domains did show a disruption in presynaptic and postsynaptic LGI1 20 signalling <sup>13</sup>. This was due to an initial decrease of Kv1 levels (at 13 days) followed by 21 AMPA receptor downregulation by 18 days post-infusion <sup>13</sup>. We also demonstrated loss of 22 Kv1.1 clusters in hippocampal CA3 (with a trend towards a reduction in glutamatergic 23 synapses) but our infusion time was shorter (7 days), and longer-term timepoints were not 24 analysed. We chose to record our EEG from the hippocampus as the expression pattern of the 25 LGI1 protein in rats shows a prominent distribution here<sup>27</sup> and LGI1 antibodies are known to 26 27 bind avidly to the CA3 area hippocampal subfield in murine sections<sup>1</sup>. Additionally, subacute reduction of hippocampal LGI1 with alternative methods (shRNA) also increases neuronal 28 network excitability<sup>28</sup>. In patients with LE and LGI-Abs, bilateral focal CA3 hippocampal 29 subfield atrophy has also been described, and the temporal lobe is the area where epileptiform 30 activity is most frequently located <sup>5, 19</sup>. The primary motor cortex has also recently been 31 identified as an additional target for LGI antibodies<sup>8</sup>. Our study did not record from motor 32 cortex however a recent elegant study using injections of dendrotoxin to block Kv1.1 33 34 channels in the rodent motor cortex faithfully recapitulated EEG changes and the clinical phenotype of tonic-dystonic seizures seen in patients with LGI1-Abs<sup>29</sup>. This alternative approach to passive transfer i.e., targeting the pathophysiological mechanisms downstream to LGI1, through a direct pharmacological inhibition of Kv1.1 channels, is a promising and exciting area in antibody-mediated seizure model development. Interestingly, although tonicdystonic seizures were not observed in the current study, some of the clinical characteristics (e.g. behavioural change, orofacial seizures) did overlap with spontaneously arising feline LGI1-autoantibody limbic encephalitis described in a recent large international cohort.<sup>30</sup>

8 In summary, in this focused study, we have developed a passive transfer *in vivo* model of 9 LGI1-mAb associated seizures which displays epileptiform activity in keeping with the 10 human disease. The use of wireless EEG telemetry and continuous video-EEG allows for 11 accurate tracking of behavioural and EEG changes. Future studies will focus on interrogating 12 new and existing data on exact pathomechanisms for *in vivo* epileptogenicity and 13 identification of possible novel epitope-specific treatment targets.

14

# **15 Data availability**

16 Data are available on request.

17

# 18 Funding

This was supported by a Wellcome Trust Fellowship [216613/Z/19/Z] to S.K.W; a senior 19 clinical fellowship from the Medical Research Council [MR/V007173/1] and Wellcome Trust 20 Fellowship [104079/Z/14/Z] to S.R.I., the German Research Foundation [FOR3004 21 SYNABS, HA6386/9-2, HA6386/10-2 to S.H. and GE2519/8-1 and GE2519/9-1 to C.G.] and 22 the European Research Council [ERC CoG 865634] to S.H., the Schilling Foundation to C.G. 23 , the German Research Foundation (SI-1969/2-1, SI-1969/3-1) and SMA Europe to C.M.S., 24 25 and by the National Institute for Health Research (NIHR) Oxford Biomedical Research Centre (BRC). For the purpose of Open Access, the author has applied a CC BY public 26 copyright licence to any Author Accepted Manuscript (AAM) version arising from this 27 submission. The views expressed are those of the author(s) and not necessarily those of the 28 29 NHS, the NIHR or the Department of Health.

30

# 1 Competing interests

SRI has received honoraria/research support from UCB, Immunovant, MedImmun, Roche,
Janssen, Cerebral therapeutics, ADC therapeutics, Brain, CSL Behring, and ONO Pharma and

4 receives licensed royalties on patent application WO/2010/046716 entitled 'Neurological

5 Autoimmune Disorders', and has filed two other patents entitled "Diagnostic method and

6 therapy" (WO2019211633 and US-2021-0071249-A1; PCT application WO202189788A1)

7 and "Biomarkers" (PCT/GB2022/050614 and WO202189788A1). SW has received honoraria

8 from UCB.

9

# 10 Supplementary material

11 Supplementary material is available at *Brain* online.

12

# 13 **References**

14

 Irani, S.R., et al., Antibodies to Kv1 potassium channel-complex proteins leucinerich, glioma inactivated 1 protein and contactin-associated protein-2 in limbic encephalitis, Morvan's syndrome and acquired neuromyotonia. *Brain*, 2010. **133**(9): p. 2734-48.

18 2. Aurangzeb, S., et al., LGII-antibody encephalitis is characterised by frequent,
19 multifocal clinical and subclinical seizures. *Seizure*, 2017. 50: p. 14-17.

Irani, S.R., et al., Faciobrachial dystonic seizures precede Lgi1 antibody limbic
encephalitis. *Ann Neurol*, 2011. 69(5): p. 892-900.

4. Thompson, J., et al., The importance of early immunotherapy in patients with
faciobrachial dystonic seizures. *Brain*, 2018. 141(2): p. 348-356.

5. Miller, T.D., et al., Focal CA3 hippocampal subfield atrophy following LGI1 VGKCcomplex antibody limbic encephalitis. *Brain*, 2017. 140(5): p. 1212-1219.

26 6. Li, L.H., et al., Clinical and electrographic characteristics of seizures in LGI127 antibody encephalitis. *Epilepsy Behav*, 2018. 88: p. 277-282.

28 7. Steriade, C., et al., Electroclinical features of seizures associated with autoimmune
29 encephalitis. *Seizure*, 2018. 60: p. 198-204.

30 8. Navarro, V., et al., Motor cortex and hippocampus are the two main cortical targets in
31 LGI1-antibody encephalitis. *Brain*, 2016. 139(Pt 4): p. 1079-93.

32 9. Steriade, C., et al., Subclinical temporal EEG seizure pattern in LGI1-antibody33 mediated encephalitis. *Epilepsia*, 2016. 57(8): p. e155-60.

Kornau, H.C., et al., Human Cerebrospinal Fluid Monoclonal LGI1 Autoantibodies
Increase Neuronal Excitability. *Ann Neurol*, 2020. 87(3): p. 405-418.

36 11. Ohkawa, T., et al., Autoantibodies to epilepsy-related LGI1 in limbic encephalitis

37 neutralize LGI1-ADAM22 interaction and reduce synaptic AMPA receptors. J Neurosci,

38 2013. **33**(46): p. 18161-74.

- Baudin, P., et al., In Vivo Injection of Anti-LGI1 Antibodies into the Rodent M1
   Cortex and Hippocampus Is Ineffective in Inducing Seizures. *eNeuro*, 2023. 10(3).
- 3 13. Petit-Pedrol, M., et al., LGI1 antibodies alter Kv1.1 and AMPA receptors changing
  4 synaptic excitability, plasticity and memory. *Brain*, 2018. 141(11): p. 3144-3159.
- 14. Ramberger, M., et al., Distinctive binding properties of human monoclonal LGI1
  autoantibodies determine pathogenic mechanisms. *Brain*, 2020. 143(6): p. 1731-1745.
- 7 15. Wright, S.K., et al., Multimodal electrophysiological analyses reveal that reduced
  8 synaptic excitatory neurotransmission underlies seizures in a model of NMDAR antibody9 mediated encephalitis. *Commun Biol*, 2021. 4(1): p. 1106.
- 16. Kreye, J., et al., Encephalitis patient-derived monoclonal GABAA receptor antibodies
  cause epileptic seizures. *J Exp Med*, 2021. 218(11).
- 12 17. Wright, S., et al., Epileptogenic effects of NMDAR antibodies in a passive transfer
  13 mouse model. *Brain*, 2015. 138(Pt 11): p. 3159-67.
- 14 18. Modebadze, T., et al., A Low Mortality, High Morbidity Reduced Intensity Status
  15 Epilepticus (RISE) Model of Epilepsy and Epileptogenesis in the Rat. *PLoS One*, 2016.
  16 11(2): p. e0147265.
- 17 19. Roberto KT, et al. Electroencephalographic findings in antileucine-rich glioma18 inactivated 1 (LGI1) autoimmune encephalitis: A systematic review. *Epilepsy Behav.* 2020
  19 112:107462
- 20. Rose, N.R., et al., Defining criteria for autoimmune diseases (Witebsky's postulates
  21 revisited). *Immunol Today*, 1993. 14 (426–430.
- 22 21. Jagtap SA, Aurangabadkar K, Joshi A, Chitnis S, Rathod M, Khade H. LGI1
  23 Encephalitis: Autoimmune Epilepsy or Movement Disorder. Ann Indian Acad Neurol.
  24 2023;26(5):831-834.
- 25 23. Fukata, Y., et al., Epilepsy-related ligand/receptor complex LGI1 and ADAM22
  26 regulate synaptic transmission. *Science*, 2006. **313**(5794): p. 1792-5.
- 27 23. Sell, J., et al., Comparative Effects of Domain-Specific Human Monoclonal
  28 Antibodies Against LGI1 on Neuronal Excitability. *Neurol Neuroimmunol Neuroinflamm*,
  29 2023. 10(3).
- Rosch, R.E., et al., NMDA-receptor antibodies alter cortical microcircuit dynamics. *Proc Natl Acad Sci U S A*, 2018. 115(42): p. E9916-e9925.
- van Sonderen, A., et al., Anti-LGI1 encephalitis: Clinical syndrome and long-term
  follow-up. *Neurology*, 2016. 87(14): p. 1449-1456.
- 26. Loscher, W., S. Cramer, and U. Ebert, Differences in kindling development in seven
  outbred and inbred rat strains. *Exp Neurol*, 1998. 154(2): p. 551-9.
- 36 27. Schulte U, Thumfart JO, Klöcker N, Sailer CA, Bildl W, Biniossek M, et al. The 37 epilepsy-linked Lgi1 protein assembles into presynaptic Kv1 channels and inhibits 38 inactivation by Kv $\beta$ 1. *Neuron*. 2006; 49: 697–706
- 28. Lugarà E, Kaushik R, Leite M, et al. LGI1 downregulation increases neuronal circuit
  40 excitability. Epilepsia. 2020;61(12):2836-2846.
- 41 29. Baudin P, Whitmarsh S, Cousyn L, et al. Kv1.1 channels inhibition in the rat motor
  42 cortex recapitulates seizures associated with anti-LGI1 encephalitis. Prog Neurobiol.
  43 2022;213:102262.
- 30. Binks S, Crawford AH, Woodhall M, et al., 067 One Health: Clinical characteristics
  of spontaneously-arising feline LGI1- autoantibody limbic encephalitis in a large
  international cohort Journal of Neurology, Neurosurgery & Psychiatry 2022;93:e2.15
- 47
- 48 49
- 50

# 1 Figure legends

Figure 1 Patient derived monoclonal LGI1 antibodies bind most strongly to CA3 region 2 3 of the hippocampus and cause a decrease of total synaptic clusters of Kv1.1 following 7day intracerebroventricular infusion. (A). Experimental timeline. (B) (i-ix) Representative 4 5 confocal images of hippocampus (CA1 (i-iii), CA3 (iv-vi) and DG (vii-ix) from sagittal brain slice prepared after 7 days of infusion of control-Ab (n=6; i, iv, and vii)) and LGI1-mAbs 6 7 (directed against epitopes EPTP n=3; (ii, v, and viii) and LRR n=4 (iii, vi, and ix)) shows typical staining pattern with secondary anti-human IgG (green). Scale bar =  $250\mu$ m. (C) (i-vi) 8 9 Example images of STED microscopy of a section from the CA3 region stained for Kv1.1 (left) and VGlut1 (middle) as well as the merge (right) from brains treated with control-Ab 10 (top) or LGI1-mAbs (bottom). (D) LGI1 monoclonal antibody treated brain slices (17 slices 11 from 7 animals; n=5 images of CA3, CA1, DG regions from each slice) median fluorescent 12 intensity log EC50 values compared to control-Ab (21 brain slices from 6 animals; n=5 13 images of CA3, CA1, DG regions from each slice ); \*P < 0.05, Mann-Whitney). (E) 14 Quantification of density of total Kv1.1 clusters in pooled analysis of CA3 region in animals 15 treated with control or LGI monoclonal antibodies. (F) Quantification of density of total 16 vGlut1 clusters in pooled analysis of CA3 region in animals treated with control or LGI 17 18 monoclonal antibodies. (G) Quantification of total KV1.1 clusters colocalizing with vGlut1. **D-G** (n= 5 images of CA3 region for each brain sample from 6 control-Ab infused and 7 19 20 LGI1-mAb infused rodents (4 LRR specific and 3 EPTP specific). Data are presented as mean  $\pm$  SEM. \**P* < 0.05, \*\**P* < 0.01, \*\*\**P* < 0.001. 21

22

Figure 2 Intracerebroventricular infusion (7 days) of patient derived monoclonal LGI1 23 antibodies induces seizures and epileptiform activity in vivo (A) Example EEG recording 24 25 from rodent using wireless EEG transmitter during intracerebroventricular infusion of LGImAb directed towards LRR (Suppl Video1) (i), highlighted/hatched EEG area expanded in 26 27 (ii); control-Ab infused rodent EEG shown for comparison in (iii). (B) Total number of 1-28 second ictal events recorded during 7-day intracerebroventricular infusion of LGI-mAbs 29 compared (n=8; 4 LRR epitope and 4 EPTP) to control-Ab infused animals (n=6). Squares 30 represent rodents infused with LRR epitope specific LGI1-mAb, circles represent EPTP 31 epitope specific LGI1-mAb. Data are presented as mean  $\pm$  SEM. (Mann-Whitney, \*P < 0.05). (C) Averaged hourly EEG coastline length for LGI1-mAb infused rodents (n=8; 4 LRR 32 epitope and 4 EPTP) compared to controls (n=6) (Mann-Whitney, \*\*\*P < 0.001). (D) Hourly 33

1 EEG power averages over 7-day infusion and recording period of LGI1-mAbs versus controls (Mann-Whitney, \*\*\*P < 0.001); rodent numbers as in (B,C). (E) Plots of frequency, duration, 2 3 and recurrence of ictal events in each animal infused with LGI1-mAbs directed towards the LRR epitope. Video-EEG analysis revealed tonic-clonic seizures (Supp V2) (i,ii,iii), non-4 5 convulsive ictal events (i,ii,iii), myoclonic jerks and wet-dog shakes (iv; Supp V3). (F) Plots of frequency, duration and recurrence of ictal events in each animal infused with LGI1-mAbs 6 directed towards the EPTP epitope. Video-EEG analysis revealed tonic-clonic seizures (i), 7 8 non-convulsive ictal events (i-iv), and myoclonic jerks (iii).

9

## 10 Figure 3 Intracerebroventricular infusion (7 days) of patient derived monoclonal LGI1

antibodies causes enhanced startle response and aggression but has no effect on memory 11 in behavioural tests. (A) Post-seizure behavioural battery scores in rodents over 7-day 12 intracerebroventricular infusion of LGI-mAbs (n=9; 5 LRR epitope and 4 EPTP) and control-13 Ab (n=6) (one-way ANOVA, \*\*P < 0.01); dotted line represents a score of 10, over this is 14 considered suggestive of behaviour consistent with spontaneous recurrent seizures. (B) Novel 15 object recognition (NOR) index calculated at three timepoints of LGI1-mAbs (n=7; 4 LRR 16 epitope and 3 EPTP) and control-Ab (n=6) infused rats. (C) Distance travelled by rats and 17 (D) velocity measured during NOR tests in rats infused with control-Ab (n=6) and LGI1-18 mAbs (n=7; 4 LRR epitope and 3 EPTP) infused rats. Data are presented as mean  $\pm$  SEM. 19

- 20
- 21 22





